

## Memo

Please share as appropriate.

To: North Dakota Clinical Laboratories

From: Christie Massen (clmassen@nd.gov)

Director, Division of Microbiology

Date: 5/14/2020

Re: FDA EUA SARS-CoV-2 IgG

The 2020 Directory of Services was sent out earlier this year. The Directory of Services is a resource meant to guide you from specimen collection through the processing and reporting of tests offered by our laboratory.

SARS-CoV-2	Acceptable Specimen: 2 ml serum or plasma	
IgG CPT CODE: 86769	<b>Stability:</b> Ambient: 16 hours; Refrigerated: 7 days; Frozen: >7 days	No Charge
	Normal Value: Negative	
	Turn Around Time: 3 days	

\*Positive results may be due to past or present infection with non-SARS-CoV-2 coronavirus strains, such as coronavirus DKU1, NL63, OC43, or 229E. Results from antibody testing should not be used as the sole basis to diagnose or exclude SARS-CoV-2 infection or to inform infection status.

Please refer to the <u>Fact Sheet for Healthcare Providers</u> and the <u>Fact Sheet for Recipients</u> for additional information.

If you have any questions, please contact our facility for the recommended guidance. **We are here to serve you!**